

The PneuTrac™ Unweighting System allows patients to perform therapeutic activities in an upright, functional position without altering gait patterns in a safe and secure environment for both client and practitioner.

A TRULY FLEXIBLE SYSTEM

- The largest range of unweighting : 0 to 300 lbs.
- 3" - 36" vertical oscillation.
- Balance and Exercise modes.
- 360 degree rotation allows multi-planar exercise.
- Patented vest with adjustability for posture and gait.
- Sturdy and simple to use analog control panel.
- Customized pricing available.
- Stand alone tracks or ceiling mount.

KEY BENEFITS

- Abilities
- Bariatric
- Orthopedic
- Strokes
- Backs
- Cardiovascular
- Respiratory ailments
- Neurological
- Safety
- Balance Training
- Gait training




Analog Control Panel

Easy to read, easy to use analog control allows precise unweighting from 0 to 300 lbs. in 1 lb. increments. Features balance and exercise modes in familiar dial format.

The system can be used over most exercise equipment, treadmills, vibration platforms and balance equipment.

All PneuWeight systems require a compressor to operate.





PneuTrac Linear



PneuTrac Trolley

Fits into most existing ceiling tracks



PneuTrac Trolley

Fits into most existing ceiling tracks

The PneuTrac XY and PneuTrac Linear are easy to assemble –NO ceiling mount necessary

**New for 2018
Pneumex PneuTrac™
With Smart Bar
Technology**

Unweighting studies on the web
www.pneumex.com

